## Aplicación

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- Palareti G. Fibrinogen/fibrin Degradation Products: Pathophysiology and Clinical Application, Fibrinolysis. 1993; 7: 60-61.
- Gaffney PJ. The Occurrence and Clinical Relevance of Fibrin Fragments in Blood, Annals New York Academy of Sciences. 1983; 407-423.
- Richmond JY, McKinney RW eds. Biosafety in Microbiological and Biomedical Laboratories, U.S. Dept. of Health and Human Services, Public Health Service, 4<sup>th</sup> Edition, 1999.

- Clinical and Laboratory Standards Institute/NCCLS. Preparation and Testing of Reagent Water in the Clinical Laboratory, Fourth Edition, CLSI/NCCLS Document C3-A4; Vol. 26 No. 22.

- Nieuwenhuizen W. A Reference Material for Harmonisation of D-Dimer Assays, Thrombosis and Haemostasis. 1997; 77 (5): 1031-1033.

- Tumanova I, Serra J. New quantitative assay for D-Dimer. Method comparison study: standardization difficulties. Communication presented at the XIII International Congress on Fibrinolysis and Thrombolysis, Barcelona, June 24-28, 1996.

- Jennings I, Woods TAL, Kirchen DP, Kitchen S, Walker ID. Laboratory D-Dimer measurement: improved agreement between methods through calibration. Thrombosis and Haemostasis. 2007; 98: 1127-1135.

- Meijer P, Haverkate K, Kluit C, de Moorfoose P, Vertruggen B, Spannagl M. A Model for the Harmonisation of Test Results of Different Quantitative D-Dimer Methods, J Thromb Haemost, 2006; 95: 567-72.

ACL™ and AcuStar™ are trademarks of Instrumentation Laboratory.  
©2008 Instrumentation Laboratory.  
Issued August 2014

## Symbols used / Verwendete Symbole / Simbolos utilizados / Symboles utilisés / Simboli impiegati / Simbolos utilizados

|                         |   |
|-------------------------|---|
| <b>CONTROL NAME</b>     | IL Control name / Namen steuern / Nombre del control / Contrôle / Controllare il nome / Nome do Controlo                              |
| <b>NAME</b>             | Name / Name / Nombre / Nom / Nome / Nome  |
| <b>TARGET DOSE</b>      | Target dose / Ziel-Dosis / Valor de referencia / Valeur médiane / Dose dell'obiettivo / Dose Alvo                                     |
| <b>ACCEPTANCE RANGE</b> | Acceptance range / Referenzbereich / Rango esperado / Limites d'acceptabilité / Intervallo di accettabilità / Intervalo de aceptación |
| <b>UNITS</b>            | Measurement units / Messeinheiten / Unidades de medida / Unités de mesure / Unità di misura / Unidades de Medida                      |